

DC1015 Scan Mode

Create mailbox name / Create folder name

- Tools
- Select 'create/change mailbox / 'create/change folder'
- Select Mailbox / folder
- Touch 'enter/change mailbox name' / 'enter/change folder name'
- Key in mailbox name (e.g Accounts)
- Save

Create password for mailbox and folder

- Tools
- Select 'create/change mailbox' / 'create/change folder'
- Select mailbox / folder
- Touch 'select mailbox' / 'select folder'
- Touch 'password'
- Modify
- Input 'new password : XXXX
- Input 'Confirmation password :XXXX
- Confirm
- Close

To Scan

- Scan
- Select desire mailbox
- Touch 'Select Mailbox
- Select feature
- Start

OPTIONAL : RETRIEVE SCANNED IMAGE

To retrieve from Email :

- Inbox

To retrieve from SMB sharing:

- My network place (win 2K/XP) or Network neighbourhood(win98/95)
- Entire network
- Microsoft windows network
- Key in 'mailbox name' (e.g Accounts)
- Key in 'password'

To retrieve from FTP:

- IE
- [FTP://001@13.199.112.162](ftp://001@13.199.112.162) (IP address)
- User name : (mailbox number)
- Password : (XXXXXX)
- Click at 'Login'
- Double click at image